

AVAPS™ and AVAPS-AE™ to VOCSN: Example Settings Conversions

Use the information on this page to aid in transitioning patient settings from AVAPS to VOCSN Volume Targeted therapy. VOCSN Volume Targeted therapy delivers a set Tidal Volume to the patient by adjusting the pressure control target. See the Breath Types and Ventilation Modes chapter of the VOCSN Clinical and Technical Manual (available at VenteLife.com/VOCSNmanual) for more information.

WARNING: These conversion settings are an example only. To ensure patient safety, verify all settings are appropriate for the patient before providing therapy.

Trilogy AVAPS		→	VOCSN	
Mode	AVAPS	→	Mode	Vol. Targeted-PS
Tidal Volume	500	→	Tidal Volume	500
IPAP Max	25	→	High Inspiratory Pressure alarm (IPAP Max+5)	30
IPAP Min	8	→	Pres. Minimum (IPAP Min-EPAP)	4
EPAP	4	→	PEEP	4
BPM	15	→	Breath Rate (or Backup Rate)	15
Ti	1.0	→	Inspiratory Time	1.0
AVAPS Rate (cmH2O per min)	1 (Slow) to 5 (Fast)	→	Pres. Adj. Rate (cmH2O per breath)	Slow (1) or Fast (3)

Trilogy AVAPS-AE		→	VOCSN	
Mode	AVAPS-AE	→	Mode	Vol. Targeted-PS
Tidal Volume	500	→	Tidal Volume	500
Maximum Pressure	30	→	High Inspiratory Pressure alarm (EPAP+Pressure Support+5)	25
Pressure Support Max	15			
Pressure Support Min	8	→	Pres. Minimum (cmH2O above PEEP/EPAP)	8
EPAP Max Pres.	8	→	PEEP	5
EPAP Min Pres.	5			
BPM	15	→	Breath Rate (or Backup Rate)	15
Ti	1.0	→	Inspiratory Time	1.0
AVAPS Rate (cmH2O per min)	0.5 (Slow) to 5 (Fast)	→	Pres. Adj. Rate (cmH2O per breath)	Slow (1) or Fast (3)

VOCSN Volume Targeted Algorithm

